



Supplementary Figure 3: Paris system intracavitary applicator. (a) Radium-226 source tubes. (b) Blind ended rubber tandem having three radium sources (1:1:0.5) i.e. 13.33 mg, 13.33 mg, 6.66 mg of radium-226 sources. (c) Metal sheath of vaginal colpostat. (d) Radium tube inside rubber tandem. (e) Two Colpostat joined contain one 13.33 mg radium-226 source in each.^[18] (f) Coronal view of whole applicator assembly inside the patient^[20]